






Characteristics and Advantages

- Catheter tube holder for temporary fixation of a low-vacuum drainage bottle
- Identical silicon material to drainage tube
- Clip for medical device with silicone material same as drainage tube
- Easy operation with simple fastening structure
- Silicone friction preventing leakage of drain pipe
- Hydro-colloid skin adhesion
- Gentle adhesion along the curved surface of the skin to reduce irritation during removal




How It Works



1 Prepare the product



Apply the adhesive tape on the affected area



Insert the drain pipe into the pipe support. Fix with the silicon band



Apply additional TPU films over the area

Features

AVAILABLE CURE FIX PRODUCTS		
CODE	SIZE	BOX QTY
CF-3 (Ø3.2)	(l) 96 × (w)41 × (h)11	50pcs × 10/box
CF-4 (Ø4.8)	(l) 96 × (w)41 × (h)11	50pcs × 10/box

SAEUM **#** MEDITEC

Tel +82 70 4278 7580 Fax +82 32 661 0805
Email marketing@saeum.net Website www.saeummeditec.com







Head Office #102-1406, 36, Bucheon-ro 198 Beon-Gil, Bucheon-si, Gyeonggi-do 14557 Rep. of Korea

Factory #102-1505, 36, Bucheon-ro 198 Beon-Gil, Bucheon-si, Gyeonggi-do 14557 Rep. of Korea